

(g)(10) Standardized API - MEDENT VERSION 23.7

**Bundle**

| Name  | Type    | Is Optional? | Description   |
|-------|---------|--------------|---|
| id    | string  | TRUE         | processId-GeneratedDateTime in UTC format           |
| meta  | instant | TRUE         | UTC format timestamp of the lastUpdate instant.     |
| type  | code    | FALSE        | searchset   |
| total | int64   | TRUE         | Total number of resources included in the response. |
| link* | url     | TRUE         | Link for the included data                          |

**Practitioner**

| Name       | Type         | Is Optional? | Description   |
|------------|--------------|--------------|---|
| identifier | Identifier   | FALSE        | Identifier will be the provider NPI if available.   |
| active     | Boolean      | TRUE         | Indicates if Provider is active                     |
| name       | HumanName    | FALSE        | Name of the Provider                                |
| telecom*   | ContactPoint | TRUE         | The Provider Home, Work, or Mobile number or Email. |
| address*   | address      | TRUE         | Provider Address                                    |

**HumanName**

| Name    | Type   | Is Optional? | Description           |
|---------|--------|--------------|-----------------------|
| use     | string | TRUE         | Value will be "usual" |
| family  | string | TRUE         | Family name           |
| given*  | string | TRUE         | Given name            |
| suffix* | string | TRUE         | Suffix                |

**ContactPoint**

| Name   | Type   | Is Optional? | Description                                    |
|--------|--------|--------------|--|
| system | string | TRUE         | Values may include "phone" or "email"          |
| value  | string | TRUE         | Phone or email                                 |
| use    | string | TRUE         | Values may include "home," "work," or "mobile" |

**Address**

| Name               | Type   | Is Optional? | Description            |
|--------------------|--------|--------------|------------------------|
| use                | string | TRUE         | Value will be "work"   |
| line*              | string | TRUE         | Address Line           |
| type               | string | TRUE         | Value will be "postal" |
| address-city       | string | TRUE         | City                   |
| address-state      | string | TRUE         | State                  |
| address-postalCode | string | TRUE         | Postal Code            |

|                 |        |      |                          |
|-----------------|--------|------|--------------------------|
| address-country | string | TRUE | Country                  |
| address-period  | Period | TRUE | Period address is in use |