

(g)(10) Standardized API - MEDENT VERSION 23.7**Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Name	Type	Is Optional?	Description
practitioner	Reference	TRUE	Provider Resource
organization	Reference	TRUE	Organization information for the PractiseID
code	CodeableConcept	TRUE	Provider Type based on urn::oid::2.16.840.1.113883.6.96
specialty	CodeableConcept	TRUE	Provider Specialty Code(s)
location	Reference	TRUE	Location Resource
healthcareService	Reference	n/a	Not Supported
telecom	ContactPoint	TRUE	Provider Contact Information
endpoint	Reference	TRUE	PractiseID Information

ContactPoint

Name	Type	Is Optional?	Description
system	string	TRUE	Values may include phone, fax, or email
value	string	TRUE	Phone or email
use	string	TRUE	Values will be work