

**(g)(10) Standardized API - MEDENT VERSION 23.5 R4****Bundle**

| <b>Name</b> | <b>Type</b> | <b>Is Optional?</b> | <b>Description</b>                                  |
|-------------|-------------|---------------------|---|
| id          | string      | TRUE                | processId-GeneratedDateTime in UTC format           |
| meta        | instant     | TRUE                | UTC format timestamp of the lastUpdate instant.     |
| type        | code        | FALSE               | searchset   |
| total       | int64       | TRUE                | Total number of resources included in the response. |
| link*       | url         | TRUE                | Link for the included data                          |

**Procedure**

| <b>Name</b> | <b>Type</b>     | <b>Is Optional?</b> | <b>Description</b>  |
|-------------|-----------------|---------------------|---|
| identifier* | Identifier      | TRUE                | UUID of the Procedure record  |
| bodySite*   | CodeableConcept | TRUE                | Body Site of the Procedure  |
| category    | CodeableConcept | TRUE                | Category description  |
| reasonCode* | CodeableConcept | TRUE                | SNOMED for the reason for procedure   |
| status      | string          | FALSE               | Values will include preparation, in-progress, not-done, on-hold, stopped, completed, or unknown |
| code        | CodeableConcept | FALSE               | Code description fo the procedure   |
| subject     | Reference       | FALSE               | References the Patient resource   |
| performed   | Period          | FALSE               | When the procedure was performed  |
| performer*  | Performer       | TRUE                | References the Practitioner resource  |

**Performer**

| <b>Name</b> | <b>Type</b> | <b>Is Optional?</b> | <b>Description</b>                    |
|-------------|-------------|---------------------|---------------------------------------|
| actor       | Reference   | FALSE               | References the Practitioner resource. |
| onBehalfOf  | Reference   | TRUE                | References the Organizaton resource.  |