

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

Procedure

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Procedure record
status	string	FALSE	Values will include preparation, in-progress, not-done, on-hold, stopped, completed, or unknown
code	CodeableConcept	FALSE	Code description fo the procedure
subject	Reference	FALSE	References the Patient resource
performed	Period	TRUE	When the procedure was performed
basedOn	Reference	TRUE	ServiceRequest Resource referenced for the Procedure