

**(g)(10) Standardized API - MEDENT VERSION 23.7****Bundle**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**Provenance**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
activity	CodeableConcept	TRUE	Values include Update, Create, Error.
agent	Agent	FALSE	Indicates who took part in an activity
recorded	instant	FALSE	When the activity was recorded
target	Reference	FALSE	References the resource the provenance file is associated with.

**Agent**

<b>Name</b>	<b>Type</b>	<b>Is Optional?</b>	<b>Description</b>
onBehalfOf	Reference	TRUE	References the Organization resource.
role*	CodeableConcept	TRUE	Values include Author or Assembler if the record was created in the EMR or if it was imported respectively
type	CodeableConcept	TRUE	Values include Primary Author, Performer, or Transmitter if the record was created in the EMR, imported into the chart, or the record was exported respectively
who	Reference	FALSE	References the Organization for records exported to a requester or items imported. This will indicate the user in MEDENT if the record originated in the EMR.