

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

DocumentReference

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the Document
status	code	FALSE	If document is closed, "completed" Otherwise the value will be "in-progress"
subject	Reference	FALSE	Coding for the problem list entry. This will be SNOMED coded.
authored	instant	TRUE	When the record id was created based on the Document Date
author*	Reference	TRUE	References the Practitioner resource
item	string	TRUE	linkID will contain the Question LOINC text will contain the Question LOINC Description Answer will be the Answer LOINC information