

(g)(10) Standardized API - MEDENT VERSION 23.7

Bundle

Name	Type	Is Optional?
id	string	TRUE
meta	instant	TRUE
type	code	FALSE
total	int64	TRUE
link*	url	TRUE

Patient

Name	Type	Is Optional?
identifier*	Identifier	TRUE
active	Boolean	FALSE
name	HumanName	TRUE
telecom*	ContactPoint	TRUE
patient	Reference	FALSE
address*	Address	TRUE

Address

Name	Type	Is Optional?
use	string	TRUE
line*	string	TRUE
type	string	TRUE
address-city	string	TRUE
address-state	string	TRUE
address-postalCode	string	TRUE
address-country	string	TRUE
address-period	Period	TRUE

HumanName

Name	Type	Is Optional?
use	string	TRUE
family	string	TRUE
given*	string	TRUE
suffix*	string	TRUE

ContactPoint

Name	Type	Is Optional?
system	string	TRUE
value	string	TRUE

use	string	TRUE
-----	--------	------

Description
processId-GeneratedDateTime in UTC format
UTC format timestamp of the lastUpdate instant.
searchset
Total number of resources included in the response.
Link for the included data

Description
UUID of the Patient Record
Based on if the Member information is set as Active or Deactivated
Name
The Patient Home, Work, or Mobile number or Email.
Patient Resource
Address information

Description
Value will be "home"
Address Line
Value will be "postal"
City
State
Postal Code
Country
Period address is in use

Description
Value will be "usual"
Family name
Given name
Suffix

Description
Value will be phone or email.
Phone or email

Value will be home or mobile