

**(g)(10) Standardized API - MEDENT VERSION 23.7****Bundle**

Name	Type	Is Optional?	Description
id	string	TRUE	processId-GeneratedDateTime in UTC format
meta	instant	TRUE	UTC format timestamp of the lastUpdate instant.
type	code	FALSE	searchset
total	int64	TRUE	Total number of resources included in the response.
link*	url	TRUE	Link for the included data

**CarePlan**

Name	Type	Is Optional?	Description
identifier*	Identifier	TRUE	UUID of the patient record
status	code	FALSE	DI/Labs/General Orders If Errored - entered-in-error If Declined - revoked If RcvdCp,RvwdCp,Notified,Closed (with an Acquired Date) - completed Else - active  Referrals If errored: entered-in-error If completed or closed: completed Else: active
intent	code	FALSE	Value will be order
category	CodeableConcept	TRUE	Value will be Evaluation procedure (procedure) - SNOMED 386053000
code	CodeableConcept	FALSE	LOINC of Order or the SNOMED for the Referral Type
subject	Reference	FALSE	References the Patient Resource.
occurrence	Period	TRUE	Scheduled Date on the order
authoredOn	DateTime	TRUE	Ordered Date
requester	Reference	TRUE	Practitioner Resource
reasonCode	CodeableConcept	TRUE	Value from the Dx associated with the Order or the Reason for Referral Text

This resource will read from Orders (Lab/Diagnostic Imaging, and General Orders) and Referrals on the p:

atient chart.